

# VARILUX® LENSES AVAILABILITY

for Eyecare Professionals

Effective **FEBRUARY 2015**



VARILUX S 4D™		LENS ID	VARILUX S FIT™		LENS ID	VARILUX S DESIGN™		LENS ID
MATERIAL	PRODUCT	INDEX	COLOR	POWER RANGE	ADD POWER*	MAX CYLINDER		
PLASTIC	1.50 Index Clear	1.50		-10.00 to +6.00	0.75 to 4.00D	-6.00		
	1.50 Index Essilor Polarized	1.50	Gray, Brown	-10.00 to +6.00	0.75 to 4.00D	-6.00		
	1.50 Index Xperio UV™ Polarized	1.50	Gray, Brown	-10.00 to +6.00	0.75 to 4.00D	-6.00		
	1.50 Index Transitions® Signature™ VII	1.50	Gray, Brown	-10.00 to +6.00	0.75 to 4.00D	-6.00		
AIRWEAR® (POLYCARBONATE)	Clear	1.59		-10.00 to +6.00	0.75 to 4.00D	-6.00		
	Essilor Polarized	1.59	Gray, Brown	-10.00 to +6.00	0.75 to 4.00D	-6.00		
	Xperio UV™ Polarized	1.59	Gray, Brown	-10.00 to +6.00	0.75 to 4.00D	-6.00		
	Transitions Signature VII	1.59	Graphite Green, Gray, Brown	-10.00 to +6.00	0.75 to 4.00D	-6.00		
THIN&LITE®	1.67 Clear	1.67		-12.00 to +9.00	0.75 to 4.00D	-6.00		
	1.67 Transitions Signature VII	1.67	Gray, Brown	-12.00 to +9.00	0.75 to 4.00D	-6.00		
	1.74 Clear	1.74		-14.00 to +13.00	0.75 to 4.00D	-6.00		
	1.74 Transitions Signature VII	1.74	Gray, Brown	-14.00 to +13.00	0.75 to 4.00D	-6.00		

NOTES: MFH=14mm

\*ADD Power in steps of 0.25D. At fitting height less than 17mm, ADD power is up to +3.50.

Thin&Lite 1.74 lenses come systematic with choice of Crizal® Prevensia™, Crizal Sapphire UV™, Crizal Avancé UV™ or Crizal Alizé UV™

All Varilux digital lenses require a hard coat. Essilor AR stacks are required when No-Glare coating is applied (for example, Crizal® No-Glare lenses).

VARILUX® PHYSIO® W3+ EYECODE™		TEMPORAL	NASAL	VARILUX® PHYSIO® W3+ FIT		TEMPORAL	NASAL	VARILUX® PHYSIO® W3+		TEMPORAL	NASAL
MATERIAL	PRODUCT	INDEX	COLOR	POWER RANGE	BASE CURVES	ADD POWER*	CYLINDER				
PLASTIC	1.50 Index Clear	1.50		-10.00 to +6.00	2.00, 3.25, 4.50, 5.75, 7.25	0.75 to 3.50, 4.00D	-6.00				
	1.50 Index Essilor Polarized	1.50	Gray, Brown	-10.00 to +6.00	2.00, 3.25, 4.50, 5.75, 7.25	0.75 to 3.50D	-6.00				
	1.50 Index Xperio UV™ Polarized	1.50	Gray, Brown	-10.00 to +6.00	2.00, 3.25, 4.50, 5.75, 7.25	0.75 to 3.50D	-6.00				
	1.50 Index Transitions® Signature™ VII	1.50	Gray, Brown	-10.00 to +6.00	2.00, 3.25, 4.50, 5.75, 7.25	0.75 to 3.50, 4.00D	-6.00				
TREXA™ (TRIVEX®)	Clear	1.53		-10.00 to +6.00	2.00, 3.25, 4.50, 5.75, 7.25	0.75 to 3.50, 4.00D	-6.00				
	Transitions Signature VII	1.53	Gray	-10.00 to +6.00	2.00, 3.25, 4.50, 5.75, 7.25	0.75 to 3.50, 4.00D	-6.00				
AIRWEAR® (POLYCARBONATE)	Clear	1.59		-10.00 to +6.00	1.75, 3.00, 4.00, 5.50, 7.25	0.75 to 3.50, 4.00D	-6.00				
	Essilor Polarized	1.59	Gray, Brown	-10.00 to +6.00	1.75, 3.00, 4.00, 5.50, 7.25	0.75 to 3.50	-6.00				
	Xperio UV Polarized	1.59	Gray, Brown	-10.00 to +6.00	1.75, 3.00, 4.00, 5.50, 7.25	0.75 to 3.50	-6.00				
	Transitions Signature VII	1.59	Graphite Green, Gray, Brown	-10.00 to +6.00	1.75, 3.00, 4.00, 5.50, 7.25	0.75 to 3.50, 4.00D	-6.00				
THIN&LITE®	Transitions® XTRActive™	1.59	Gray	-10.00 to +6.00	1.75, 3.00, 4.00, 5.50, 7.25	0.75 to 3.50, 4.00D	-6.00				
	1.60 Clear	1.60		-12.00 to +8.00	1.75, 3.00, 4.00, 5.50, 7.00, 8.00	0.75 to 3.50, 4.00D	-6.00				
	1.60 Transitions Signature VII	1.60	Gray	-12.00 to +8.00	1.75, 3.00, 4.00, 5.50, 7.00, 8.00	0.75 to 3.50, 4.00D	-6.00				
	1.67 Clear	1.67		-14.00 to +9.00	1.75, 2.75, 3.75, 5.25, 6.50, 8.00	0.75 to 3.50, 4.00D	-6.00				
	1.67 Transitions Signature VII	1.67	Gray, Brown	-14.00 to +9.00	1.75, 2.75, 3.75, 5.25, 6.50, 8.00	0.75 to 3.50, 4.00D	-6.00				
	1.74 Clear	1.74		-18.00 to +13.00	1.75, 2.75, 3.75, 5.25, 6.50, 8.00, 10.00	0.75 to 3.50, 4.00D	-6.00				
1.74 Transitions Signature VII	1.74	Gray	-18.00 to +13.00	1.75, 2.75, 3.75, 5.25, 6.50, 8.00, 10.00	0.75 to 3.50, 4.00D	-6.00					

NOTES: MFH=14mm

\*ADD Power in steps of 0.25D. At fitting height less than 17mm, ADD power is up to +3.00.

Thin&Lite 1.74 lenses come systematic with choice of Crizal® Prevensia™, Crizal Sapphire UV™, Crizal Avancé UV™ or Crizal Alizé UV™

All Varilux digital lenses require a hard coat. Essilor AR stacks are required when No-Glare coating is applied (for example, Crizal® No-Glare lenses).

# VARILUX® LENSES AVAILABILITY

for Eyecare Professionals

Effective **FEBRUARY 2015**



VARILUX® COMFORT® W2+ EYECODE™			VARILUX® COMFORT® W2+ FIT			VARILUX® COMFORT® W2+			
			TEMPORAL	NASAL				TEMPORAL	NASAL
			E	W				W	W
MATERIAL	PRODUCT	INDEX	COLOR	POWER RANGE	BASE CURVES	ADD POWER*	CYLINDER		
PLASTIC	1.50 Index Clear	1.50		-10.00 to +6.00	2.00, 3.25, 4.50, 5.75, 7.25	0.75 to 3.50D	-6.00		
	1.50 Transitions® Signature™ VII	1.50	Gray, Brown	-10.00 to +6.00	2.00, 3.25, 4.50, 5.75, 7.25	0.75 to 3.50D	-6.00		
	1.50 Transitions® XTRActive®	1.50	Gray	-10.00 to +6.00	2.00, 3.25, 4.50, 5.75, 7.25	0.75 to 3.50D	-6.00		
AIRWEAR® (POLYCARBONATE)	Clear	1.59		-10.00 to +6.00	1.75, 3.00, 4.00, 5.50, 7.00	0.75 to 3.50D	-6.00		
	Essilor Polarized	1.59	Gray, Brown	-10.00 to +6.00	1.75, 3.00, 4.00, 5.50, 7.00	0.75 to 3.50D	-6.00		
	Xperio UV™ Polarized	1.59	Gray, Brown	-10.00 to +6.00	1.75, 3.00, 4.00, 5.50, 7.00	0.75 to 3.50D	-6.00		
	Transitions Signature VII	1.59	Graphite Green, Gray, Brown	-10.00 to +6.00	1.75, 3.00, 4.00, 5.50, 7.00	0.75 to 3.50D	-6.00		
	Transitions XTRActive	1.59	Gray	-10.00 to +6.00	1.75, 3.00, 4.00, 5.50, 7.00	0.75 to 3.50D	-6.00		
THIN&LITE®	1.60 Clear	1.60		-12.00 to +8.00	1.75, 3.00, 4.00, 5.50, 7.00, 8.00	0.75 to 3.50D	-6.00		
	1.60 Transitions Signature VII	1.60	Gray	-12.00 to +8.00	1.75, 3.00, 4.00, 5.50, 7.00, 8.00	0.75 to 3.50D	-6.00		
	1.67 Clear	1.67		-14.00 to +9.00	1.75, 2.75, 3.75, 5.25, 6.50, 8.00	0.75 to 3.50D	-6.00		
	1.67 Transitions Signature VII	1.67	Gray, Brown	-14.00 to +9.00	1.75, 2.75, 3.75, 5.25, 6.50, 8.00	0.75 to 3.50D	-6.00		

NOTES: MFH=14mm

\*ADD Power in steps of 0.25D. At fitting height less than 17mm, ADD power is up to +3.00.

All Varilux digital lenses require a hard coat. Essilor AR stacks are required when No-Glare coating is applied (for example, Crizal® No-Glare lenses).


## REQUIRED MEASUREMENTS FOR DISPENSING

	PD and FITTING HEIGHT	PANTOSCOPIC TILT	WRAP ANGLE	VERTEX DISTANCE	EYE ROTATION CENTER DISTANCE	NATURAL HEAD POSTURE	LEADING DOMINANT EYE	H/E RATIO	STABILITY COEFFICIENT	READING DISTANCE
VARILUX S 4D*	●	●	●		●	●	●	▲	▲	▲
VARILUX S FIT	●	●	●	●						
VARILUX S Design	●									
VARILUX PHYSIO W3+ eyecode*	●	●	●		●	●				
VARILUX PHYSIO W3+ FIT	●	●	●	●						
VARILUX PHYSIO W3+	●									
VARILUX COMFORT W2+ eyecode*	●	●	●		●	●				
VARILUX COMFORT W2+ FIT	●	●	●	●						
VARILUX COMFORT W2+	●									

\* Varilux S 4D lenses, Varilux Physio W3+ eyecode lenses and Varilux Comfort W2+ eyecode lenses are exclusive to the Visioffice® System.

## ENGRAVING LOCATIONS

The illustration below shows example of engraving locations for Varilux Physio W3+ eyecode



Refer to availability chart for each product's engraving symbol